

Embedding

- 1) Put the specimen on the strainer, and immerse it in the coolant (Fig. 1).

Coolants: Hexane-dry ice (approx. -75°C)

Hexane cooled with a cooling apparatus (-94°C)

Pentane in liquid nitrogen (approx. -130)

Isopentane in liquid nitrogen (approx. -160°C)

Note: Prepare more than 500ml of hexane at least when the hexane-dry ice is used as a coolant.

- 2) Put the proper quantity of embedding medium in the stainless steel container, then partly freeze the embedding medium (Fig. 2).

Note: Use the cooled embedding medium kept in a refrigerator.

- 3) Place the frozen specimen in the embedding medium.

Note: Place the sectioning side face down.

- 4) Submerge the container into the coolant (Fig. 3A)

Note: In order to avoid melting the specimen, carry out steps 3) and 4) rapidly.

- 5) After 5-15 seconds (the time depends on the container size), raise the embedding medium so that the top is exposed from the coolant and freeze it completely (Fig. 3B).

Note: If the embedding medium is frozen completely in the coolant, cracks will occur in the frozen specimen block.

- 6) After freezing, place the frozen specimen block on the extractor (a wooden stage) and remove it with a wooden hammer (Fig. 4).

- 7) Slightly melt the specimen holder side of the frozen block and attach it to the holder (Fig. 5).

- 8) Fix the holder to the Cryomicrotome specimen stage.



Fig. 1. Freezing the specimen.




Fig. 2. Embedding




Fig. 3. Freeze-embedding.

(A): Early freezing

(B): Complete freezing.




Fig. 4. Removing the frozen block from the container.




Fig. 5.

(A): Frozen specimen block

(B): The block fixed to the specimen holder.